

Pan-ascomycete mold, mouse monoclonal antibody [clone IE3]

Product code	ID2560-0025 / ID2560-0100
Unit size	25 μ L / 100 μ L
Concentration	See batch specific Certificate of Analysis
Alternative name	No
Clone	IE3
Host	Mouse
Antigen	25 to >250kDa glycoprotein
Isotype	IgM
Purification	Affinity purification
Applications	ELISA, WB
Recommended dilution	1:1000
Optimisation	Optimal dilutions to be determined by end user
Species reactivity	All filamentous ascomycete species
Storage buffer	Phosphate buffered saline pH7.2 with 0.095% (w/v) sodium azide
Shipping	Blue ice
Storage temperature	Store as supplied at +2°C ~ +8°C for up to 1 year