

***Aspergillus*, mouse monoclonal antibody [clone EF9]**

Product code	ID2600-0025 / ID2600-0100
Unit size	25 μ L / 100 μ L
Concentration	See batch specific Certificate of Analysis
Alternative name	No
Clone	EF9
Host	Mouse
Epitope/Antigen	β 1,5-galactofuranose (Galf) on 50 to >200kDa mannoprotein
Isotype	IgG2b (<i>in vitro</i> class switch of ED1 [ID2595]: JF5 \rightarrow ED1 \rightarrow EF9)
Purification	Protein A
Applications	ELISA, WB, IF, IEM
Recommended dilution	1:1000
Optimisation	Optimal dilutions to be determined by end user
Species reactivity	<i>Aspergillus</i> and <i>Penicillium</i> species
Storage buffer	Phosphate buffered saline pH7.2 with 0.095% (w/v) sodium azide
Shipping	Blue ice
Storage temperature	Store as supplied at +2°C ~ +8°C for up to 1 year

UNLESS OTHERWISE STATED IN OUR CATALOG OR OTHER COMPANY DOCUMENTATION ACCOMPANYING THE PRODUCT(S), OUR PRODUCTS ARE INTENDED FOR **RESEARCH USE ONLY** (RUO) AND ARE NOT TO BE USED FOR ANY OTHER PURPOSE, WHICH INCLUDES, BUT IS NOT LIMITED TO, UNAUTHORISED COMMERCIAL USES, *IN VITRO* DIAGNOSTIC USES, *EX VIVO* OR *IN VIVO* THERAPEUTIC USES OR ANY TYPE OF CONSUMPTION OR APPLICATION TO HUMANS OR ANIMALS.